


Synthetic Cannabinoids-3 Assay

Immunoanalysis and Specialty Diagnostix now offer a third homogeneous enzyme immunoassay for the detection of synthetic cannabinoids in urine. While the previous Synthetic Cannabinoid Assays target the first and second generation compounds and their respective metabolites, the new Synthetic Cannabinoid-3 Assay is intended for the analysis of AB-PINACA, AB-FUBINACA and their derivatives.

Assay Characteristics

- Qualitative detection of AB-PINACA, AB-FUBINACA, and other currently relevant synthetic cannabinoid compounds
- Excellent specificity and sensitivity
- Liquid, stable, ready-to-use reagents, calibrators and controls – no reconstitution required
- Application protocols available for all major clinical chemistry analyzers
- Packaging tailored to your laboratory's specific needs



Background

AB-PINACA was first identified as a component of synthetic cannabinoid products in 2012. It is based on an indazole backbone, whereas most of the JWH compounds have an indolyl base. The biochemical, physiological, and toxicological properties of AB-PINACA and its derivatives have yet to be determined.


Originally, AB-PINACA was developed as an analgesic medication by Pfizer in 2009. In rat discrimination studies, AB-PINACA was found to act as a potent agonist for the CB₁ and CB₂ receptors and as a full substitute for Δ⁹-THC, while being 1.5 times more potent.

AB-PINACA is sold under various street names, the most common being Cloud 9,

Crown, and Relax. Sprayed on herbal blends or as a liquid in eyedropper bottles, it is usually used with vaporizing devices, such as hookahs or e-cigarettes.

The most commonly reported adverse effects of herbal blends supposedly containing AB-PINACA include difficulty in breathing, loss of motor control, altered mental status, loss of consciousness, convulsions, coma and death. There have been a number of reported cases of deaths and hospitalizations connected to AB-PINACA.

According to a recent report of the Global Drug Survey, synthetic cannabinoids are now considered to be the most dangerous of all recreational drugs.



Cross-Reactivities

Analyte	Concentration (ng/mL)	Cross-Reactivity (%)	Analyte	Concentration (ng/mL)	Cross-Reactivity (%)
AB-PINACA Pentanoic Acid	10	100	AM2201 N-(4-hydroxypentyl)	100 000	< 0.05
AB-PINACA	10	100	JWH-018	35 000	< 0.05
AB-PINACA N-(4-hydroxypentyl)	8	125	JWH-018 N-(5-hydroxypentyl)	90 000	< 0.05
AB-PINACA N-(5-hydroxypentyl)	10	100	JWH-200	40 000	< 0.05
5-fluoro AB-PINACA	8	125	JWH-201	60 000	< 0.05
5-fluoro ABICA	20	50	JWH-250	60 000	< 0.05
5-fluoro ADBICA	12	83	JWH-250 N-(5-hydroxypentyl)	25 000	< 0.05
5-fluoro AB-PINACA N-(4-hydroxypentyl)	15	67	JWH-250 N-(4-hydroxypentyl)	30 000	< 0.05
5-fluoro ADB-PINACA	9	111	JWH-250 N-(5-carboxypentyl)	100 000	< 0.05
5-chloro AB-PINACA	15	67	A-796260	3 000	0.33
ADB-PINACA	15	67	A-834735	5 000	0.20
ADB-PINACA Pentanoic Acid	7	143	AB-005	50 000	0.20
ADB-PINACA N-(4-hydroxypentyl)	8	125	AM-2232	100 000	< 0.05
ADB-PINACA N-(5-hydroxypentyl)	6	167	UR-144	8 000	0.13
AB-FUBINACA	10	100	UR-144 N-(5-bromopentyl)	8 000	0.13
ADB-FUBINACA	10	100	UR-144 N-(5-hydroxypentyl)-β-d-glucuronide	3 000	0.33
ADBICA	20	50	UR-144 N-(5-chloropentyl)	5 000	0.20
ADBICA N-Pentanoic Acid	15	67	UR-144 N-(5-hydroxypentyl)	3 500	0.29
ADBICA N-(4-hydroxypentyl)	15	67	XLR-11	3 500	0.25
ADBICA N-(5-hydroxypentyl)	12	83	XLR-11 N-(4-hydroxypentyl)	3 000	0.33
AD-CHMINACA	15	67	XLR-11 N-(4-pentyl)	8 000	0.13

Method Comparison

		LC-MS/MS	
		+	-
HEIA	+	40	0
	-	0	40

Sensitivity: 100%
Specificity: 100%
Accuracy: 100%

Ordering Information

Reagents	Size	Order No.
K2-III (Synthetic Cannabinoids-3) Assay	25 mL	350-0025EX
K2-III (Synthetic Cannabinoids-3) Assay	60 mL	350-0060WEX
K2-III (Synthetic Cannabinoids-3) Assay	100 mL	350-0100EX
Calibrators		
AB-PINACA Pentanoic Acid, 10 ng/mL	1 x 5 mL	10030
Controls		
AB-PINACA Pentanoic Acid Control Set, 5 ng/mL Low & 15 ng/mL High	2 x 5 mL	3012

Specialty Diagnostix

Specialty Diagnostix GmbH
Sailerwöhr 16
94032 Passau, Germany
Phone +49 (0)851 988 4930 0
Fax +49 (0)851 988 4930 8
info@specialtydiagnostix.de
www.specialtydiagnostix.de